

# Santai Adaptor Kit for 10 kg Flash Columns



Santai Technologies has recently launched its largest Flash column, which is pre-packed with 10 kg silica. Loading up to 1 kg of the sample to the column in just one injection is now available with this large Flash column. To connect the Flash column with most commonly used large flow rate Flash chromatography systems available on the market, there are several adaptors for different tubings with various outer diameter (OD), including 1/2-inch, 3/8-inch, 3/16-inch and 1/8-inch.



With these adaptors, the large Flash column can be connected to any large flow rate Flash chromatography systems available on the market, where tubings with different OD are used.

*When connecting with 1/2-inch OD tubing:*



*When connecting with 3/8-inch OD tubing:*



*When connecting with 3/16-inch OD tubing:*



***When connecting with 1/8 inch-OD tubing:***



\*Please be noted that the adaptor for the inlet and outlet of the 10 kg Flash column are the same so only the inlet adaptor is shown in the pictures.

#### **Ordering Information**

Item No.	Description
BLL-NPT635-127	10kg cartridge adaptor for 1/2-inch OD tubing.
BLL-NPT635-095	10kg cartridge adaptor for 3/8-inch OD tubing.
BLL-NPT635-048	10kg cartridge adaptor for 3/16-inch OD tubing.
BLL-NPT635-032	10kg cartridge adaptor for 1/8-inch OD tubing.

---

## **Santai Science Inc.**

**Add :** 214 Brunswick, Pointe-Claire, Montréal, QC H9R 1A6, Canada

**Tel :** +1 514 505 1378

**Email:** [order@santaisci.com](mailto:order@santaisci.com)

[support@santaisci.com](mailto:support@santaisci.com)